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1. Houts, G.E., Miyagi, M., Ellis, C., Beard, D., and Beard, J.W. (1979) *J.Virol.* 29(2):517-522.
2. Guide to Molecular Cloning Techniques. Methods in Enzymology, Volume 152. pp 316-325. Edited by Shelby Berger and Alan R. Kimmel. Academic Press, Inc.

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
Store at -20 °C

Protocol

1. RT reaction
 - a. Reaction setup for cDNA synthesis

5 x RT Buffer	2 μl
dNTP Mixture (10Mm)	1 μl
Oligo-dT	
or Random Hexamer Primer**	0.5 μl
or special downstream primer	
RNase inhibitor (40U/ul)	0.5 μl
AMV Reverse Transcriptase (5U/ul)	0.5 μl
RNA Control or Sample RNA*	x μl
RNase free H ₂ O	5.5-x μl
total	10 μl

Notes

*Please add 20ul RNase free H₂O before using RNA Control and dissolve sufficiently. Add 2ul/test when you do controls. Sample RNA valume can be add up to 5.5ul (≤1ug).

b、 RT reaction condition

(room temperature, 10min) **

42-60°C ***, 45min

95°C ****, 5min

ice bath, 5min

cDNA for PCR reaction.

Notes

** If using random Hexamer Primer, holding at room temperature for 10 minutes.

*** Reaction temperature can be elevated for RNA template with second structure, lower than 60°C, or holding at 70°C for 10 minutes after adding RNA template and primers, then holding on ice for 5 minutes.

**** Heat to 95°C, let AMV Reverse Transcriptase denature.

2. PCR protocol

a、 Reaction setup for PCR

10 x PCR Buffer (include 15mM Mg ²⁺)	2.5 µl
dNTP Mixture (10mM)	0.5 µl
Upstream primer (5 µM)	0.5 µl
Downstream primer (5 µM)	0.5 µl
Taq mix DNA polymerase	0.5 µl
cDNA*	2.5 µl
ddH ₂ O	18 µl
total	25 ul

b、 PCR reaction condition

94°C 3min

94°C 30s

37°C-65°C ** 30s

72°C **** 45s-7min

72°C 5min

} 25-35cycles ***

Note:

* RT products can up to 5.0 µl.

** 5°C lower than T_m value of primers; Control primer is 55°C.

*** Add to 40-45 cycles when detecting rare RNA templates .

**** usually 1kb/min . Control primer is 45s.

c、 Analyze 3 - 5µl of the reaction products by agarose gel electrophoresis.

Using Tips

1. Total RNA or mRNA can be used as RNA template, and we suggest using Biozol(BSC51M1)to isolate high quality RNA;
2. RNase contamination should be avoided and follow the measures:
 - 1) Ware one time gloves and respirator because of the RNases in saliva and skin;
 - 2) Use special instruments and consumables and handle in specific areas;
 - 3) Consumables should be treated at 180°C for 60min or 37°C for 12 hours with 0.1% DEPC H₂O followed by sterilized at 121°C for 30min;
3. AMV reverse transcriptase, Taq polymerase and RNase inhibitor should be slowly pipetted after centrifuging, and put back at -20°C later;
4. Avoid frequently freezing and thawing dNTP;
5. Specific primers should be used and concentration of primers should be optimized and We suggest 0.4μM as a starting point for optimizing; Oligo-dT and Random primers are not suitable for this kit;
6. Concentration of MgCl₂ can be optimized and increasing 0.5mM/time is suggested when the length of target template is longer than 2kb;
7. Quality of primers affect the performance of BioRT One Step RT-PCR Kit; Factors such as GC percentage, length of primers and site of primers should be considered and we suggest primers be designed using special soft wares.

Required materials for RT-PCR

Instruments and consumables	Reagents
Centrifuge	DEPC
Pipettes	ddH ₂ O
RNase free 1.5ml tubes	Electrophoresis Buffer
Water bath instruments	Loading Buffer
Gel electrophoresis instruments	DNA Marker
PCR tubes	
Tips	

Reference


1. Houts, G.E., Miyagi, M., Ellis, C., Beard, D., and Beard, J.W. (1979) *J. Virol.* 29(2):517-522.
2. Guide to Molecular Cloning Techniques. Methods in Enzymology, Volume 152. pp 316-325. Edited by Shelby Berger and Alan R. Kimmel. Academic Press, Inc.

Description

RT PCR is a useful technique for the investigation of gene expression, viral load, pathogen detection, and numerous other applications. When analyzing gene expression or viral load, the RNA of interest first needs to be reverse transcribed into cDNA. The subsequent PCR can be performed separately (two-step RT-PCR).

BioRT Two Step RT-PCR Kit is designed for two step RT-PCR of RNA samples from various sources. The kit includes all the necessary reagents for cDNA synthesis and subsequent PCR. Either total RNA, messenger RNA, viral RNA or in vitro transcribed RNA can be used as a template for reverse transcription. The kit includes both random primers and oligo(dT18) primers. The user can choose either of these or alternatively use gene specific primers. The reverse transcriptase in the BioRT Two Step RT-PCR Kit is AMV, which provides up to 60°C RT temperature and provides higher sensitivity and higher yield to cDNA's synthesis and PCR. The performance of the PCR step is based on a high fidelity Taq mix DNA polymerase. The kit can synthesis up to 12kb cDNA and 6kb DNA. The reaction buffer is optimized for a kind of reaction system (10ul for RT and 25ul for PCR).

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