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Website: www.bioer.com.cn

3. AMV reverse transcriptase, Taq polymerase and RNase inhibitor should be slowly pipetted after centrifuging. The remainder is better to store at -20°C as soon as possible;
4. Avoid frequently freezing and thawing dNTP;
5. Specific primers should be used and concentration of primers should be optimized and We suggest 0.4µM as a starting point for optimizing; Oligo-dT and Random primers are not suitable for this kit;
6. 2mM MgCl₂ is already included in 10*RT-PCR buffer. The concentration of MgCl₂ can be adjusted depending on different conditions. 2.5-3mM is suggested when the length of target template is longer than 2kb;
7. The primer design directly effects the performance of RT-PCR. Factors such as GC percentage, length of primers and site of primers sites should be considered. We suggest using software, such as Primer Premier5.0.

Protocol

Prepare the reaction mix by combining the indicated components

10 x RT-PCR Buffer	2.5 µl
dNTP Mixture	4 µl
RNase inhibitor	1 µl
Sense primer(10 µM)	1 µl
Antisense primer(10 µM)	1 µl
AMV reverse transcriptase	0.5 µl
Taq polymerase	0.5 µl
Total RNA <1ug	X µl
RNase free H ₂ O	14.5-X µl
Total Volume	25 µl

The reaction system can be enlarged to 50µl in accordance with the corresponding proportion.

REACTION CONDITIONS

42°C*	30min	} 30-40cycles***
94°C	5min	
94°C	30s	
50-65°C**	30s	
72°C	1kb/min	
72°C	5min	

Analyze 5µl of the reaction products by agarose gel electrophoresis

Notes

- * Reaction temperature can be elevated moderately for RNA template with second structure
- ** 5°C lower than T_m value of primers
- *** Add to 45-50 cycles when detecting rare RNA templates

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
Store at **-20 °C**

Required materials for RT-PCR

Instruments and consumables	Reagents
Centrifuge	DEPC
Pipettes	ddH ₂ O
RNase free 1.5ml tubes	Electrophoresis Buffer
Water bath instruments	Loading Buffer
Gel electrophoresis instruments	DNA Marker
PCR tubes	
Tips	

Usage Tips

- Total RNA or mRNA can be used as RNA template, and we suggest using Biozol(BSC51M1)to isolate high quality RNA;
- RNase contamination should be avoided:
 - 1) Wear one-off gloves and respirator because of the Rnases in saliva and skin;
 - 2) Use special instruments and consumables and handle in specific areas;
 - 3) Consumables should be dry heat sterilized at 180°C for 60min, or treat at 37°C with 0.1% DEPC aqueous solution for 12 hours followed by sterilization at 121°C for 30min.


Description

RT-PCR (Reverse Transcription-Polymerase Chain Reaction) is a technique in which an RNA strand is "reverses" transcribed into its DNA complement, followed by amplification of the resulting DNA using a polymerase chain reaction (PCR). RT-PCR can be used to examine gene expression level in cells and tissues, clone the specific gene of cDNA sequences and test RNA viruses.

Unlike two steps RT-PCR methods, in which RT and PCR separated in two tubes, BioRT One Step RT-PCR Kit is designed to combine RT and PCR in the same response system. No additional reagents are needed during the process. Our kits adopt High-quality AMV reverse transcriptase produced in USA and Taq polymerase. We introduce an optimized buffer to ensure the maximize effectiveness of these two enzymes. Cooperation of AMV reverse transcriptase and highly purified recombinant Taq polymerase can obtain 4kb extent of products

BioRT One Step RT-PCR Kit adopts one tube system. Because operator doesn't need to open the lid during the reaction process, this user-friendly improved version avoids cross contamination and has the merits of high specificity and sensitivity.

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