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Important Parameters

Template

Using $10\text{-}10^7$ copies of plasmid DNA, 100pg to 1ug of genomic DNA or 1pg to 100ng of cDNA. In order to acquire the best experiment result, the amplification length should be approximately 80-500bp.

Primers

Primer design is one of the most important parameters when using SYBR Green detection system. We strongly recommend using primer design software such as Oligo, Primer5. When designing primers, the primers length should be approximately 15-40bp. Optimal results may require a titration of primer concentrations between 100 and 500 nM. A final concentration of 200 nM per primer is effective for most reactions.

Magnesium Concentration

MgCl₂ is included in the SYBR Mix at a final concentration of 2 mM. This concentration works well for most targets; however, the optimal concentration may range from 2 to 6 mM. If necessary, use the 25 mM MgCl₂ provided in the kit to increase the concentration, as shown below.

Final MgCl ₂ concentration	Add this volume of MgCl ₂
3.0 mM	2 ul
4.0 mM	4 ul
5.0 mM	6 ul
6.0 mM	8 ul

ROX Reference Dye

Rox Reference Dye can be included in the reaction to normalize the fluorescent reporter signal, for instruments that are compatible with that option. Use the following table to determine the amount of Rox Reference Dye to use with a particular instrument:

Intrument	Amount of ROX Reference Dye (50×) per 50 ul reaction (ul)
ABI7000, 7300, 7300, 7900HT, 7900HT Fast	1
ABI7500, Stratagene Mx3000, Mx3005P, and Mx4000	0.1

Using of the kit

- 1 Defrosting the reagents (2×SYBR Mix、ddH₂O) in the kit at room temperature. Before preparing PCR reagents, please centrifuge all reagents for a few seconds.
- 2 Set up reactions as specified below. Volumes for a single 50ul reaction are listed.

Component	Single (ul)
2×SYBR Mix(with 4.0mM Mg ²⁺)①	25
PCR Forward Primer (10uM)	1
PCR Reverse Primer (10uM)	1
Taq Polymerase	0.3
ddH ₂ O	18.7
Template②	4
Total	50

① 2×SYBR Mix includes PCR buffer, Mg²⁺, dNTP mixture, SYBR Green I etc.

② According to different copies of target Gene between different kinds of DNA template, if possible user should confirm the additional quantity of DNA template by doing Gradient dilution.

3 After mixing PCR reagents above, distribute them into 0.2ml PCR tubes. Add the samples as per 4-10ul each into PCR tubes above. .

4 PCR Reaction

Program your Real Time PCR as shown below. Optimal temperatures and incubation times may vary.

<p>Three step standard cycling program</p> <p>94°C, 2min ;</p> <p>94°C, 10s ;</p> <p>60°C, 15s ;</p> <p>72°C, 30s, 40cycles</p> <p>Choose F1 channels when collecting fluorescent signals and set fluorescent signals detecting at 72°C</p>	<p>Two step standard cycling program</p> <p>94°C, 2min ;</p> <p>94°C, 10s ;</p> <p>60°C, 30s; 40cycles</p> <p>Choose F1 channels when collecting fluorescent signals and set fluorescent signals detecting at 60°C</p>
---	--

Before running, please adjust gain to make the F1 (SYBR Green I) background ≤20.

Program your melting curve as shown below.

- 95°C, 2min;
- 72°C, 1min;
- 95°C, 30s, step 0.5°C/s;
- 30°C, 1min.

Note: when user is selecting the PCR programme, the three step methods will be recommended for amplification. If the amplification result of two-step is the same as the one of three-step method, user also can use two-step method instead of three-step method for saving time.

Result analysis and judgments

1 Input the concentrations of four positive References. Select fit point method to analyze. Confirm the base line (zero adjustment) by getting the fluorescent signals.

Make the noise limit just beyond the peak of the amplification curve (rule less noise line) of normal negative control; then do quantitative analysis. You could also adjust by yourself according to the condition of instrument's noises.

- 2 Please record the concentration or CT value of non-known sample after analysis.
- 3 There are auto and manual for users choosing in melting curve analysis and record T_m value after analysis.
- 4 Users also can use agarose gel electrophoresis to analyze specificity of PCR products

Note:

1. Please read this manual carefully before beginning the experiment.
2. Don't touch hand directly when you use one-time consumable such as pipette and centrifuge tube.
3. Be caution during the whole process to ensure the accuracy of the experiment. After finishing the experiment, in order to avoid pollution, please use 75% alcohol to clear the worktable.
- 4 The melted reagent should be stored as shorten as possible, and store them in environment of -20°C immediately after testing.
5. Please store reagent or PCR reaction system in shadow place due to the fluorescence dye contained
6. Nucleic acid is easy to be contaminated and resulted in to the trusted result, so user should strictly do the amplification experiment according to standard PCR process.

Sample storage

Please store the sample at -20°C and avoid repeating frozen.

Sample transportation

Please transport the sample at or lower than 4°C.

TECHNICAL SUPPORT:

For technical support, please dial phone number +86-571-87774567-5211 or 87774575, by fax to +86-571-87774303, or by email to reagent@bioer.com.cn.

CONTACT INFORMATION OF MANUFACTURER.

HANGZHOU BIOER TECHNOLOGY CO.,LTD.

Address: 1192 Bin An Rd., Hi-Tech (Binjiang) District, Hangzhou, 310053, China

Website: www.bioer.com.cn.

TEL: +86-571-87774575;

FAX:+86-571-87774565.

Description

Adopting nucleic acid amplification technique and combining intermediate fluorescence method, SYBR Green I real-time PCR can do detection and quantification on target gene fast and properly. The kit offers 2×SYBR Mix, including all components except for Taq DNA polymerase and primers. During operation, user can do real-time PCR reaction directly by adding primer, template and Taq DNA polymerase, which is very simple.

SYBR Green I is one kind of fluorescence dye, which will excited fluorescence signal after specially inserted in dsDNA molecule. In PCR reaction system, when SYBR Green I dye combine with dsDNA, and excited with excited light, via testing strength of fluorescence, user can calculate the quantity of dsDNA. The advantage of fluorescence dye is that: it can be used in testing any dsDNA amplification, without design of probe, and make the testing method easier, meanwhile, it can decrease cost, what's more, user can analyze melt curve directly. As soon as the PCR test finish. By the analysis of melt curve, user can judge whether there are existing aberrance or non-specific amplification. The Max Absorb wave length of SYBR Green I is 497nm, the Max. Emission wavelength is about 520nm.

Component

component	size/ul	Quantity
2×SYBR Mix(with 4.0mM Mg ²⁺)	1250	4
Taq DNA Polymerase	32	2
25mM MgCl ₂	1000	2
ddH ₂ O	1000	4
ROX Reference (50×)	200	1

Instrument

Bioer Line-Gene Real Time PCR Detection System, or the same kind of the instruments in other companies.

Quality Control

This product is tested functionally in Real Time using plasmid DNA. Kinetic analysis must demonstrate a linear dose response with decreasing target concentration over four orders of magnitude, and the relative parameter of standard curve should be ≤ -0.980 .

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