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1. Houts, G.E., Miyagi, M., Ellis, C., Beard, D., and Beard, J.W. (1979) *J.Virol.* 29(2):517-522.
2. Guide to Molecular Cloning Techniques. Methods in Enzymology, Volume 152. pp 316-325. Edited by Shelby Berger and Alan R. Kimmel. Academic Press, Inc.

BioRT cDNA First Strand Synthesis Kit

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
Store at -20 °C

Protocol
RT reaction

a、Reaction setup for cDNA synthesis

5 x RT Buffer	2 μl
dNTP Mixture (10mM)	1 μl
Oligo-dT	
or Random Hexamer Primer**	0.5 μl
or special downstream primer	
RNase inhibitor (40U/ul)	0.5 μl
AMV Reverse Transcriptase (5U/ul)	0.5 μl
Sample RNA*	x μl
RNase free H ₂ O	5.5-x μl
total	10 μl

Notes

* When you do controls. Sample RNA valume can be add up to 5.5ul (≤1ug).

b、RT reaction condition

(room temperature, 10min) **

42-60°C ***, 45min

95°C ****, 5min

ice bath, 5min

cDNA for PCR reaction.

Notes

** If using random Hexamer Primer, holding at room temperature for 10 minutes.

*** Reaction temperature can be elevated for RNA template with second structure, lower than 60°C.

**** Heat to 95°C let AMV Reverse Transcriptase denature.

Using Tips

1. Total RNA or mRNA can be used as RNA template, and we suggest using Biozol(BSC51M1)to isolate high quality RNA;
2. RNase contamination should be avoided and follow the measures:
 - 1) Ware one time gloves and respirator because of the RNases in saliva and skin;
 - 2) Use special instruments and consumables and handle in specific areas;
 - 3) Consumables should be treated at 180°C for 60min or 37°C for 12 hours with 0.1% DEPC H₂O followed by sterilized at 121°C for 30min;
3. AMV reverse transcriptase, Taq polymerase and RNase inhibitor should be slowly pipetted after centrifuging, and put back at -20°C later;
4. Avoid frequently freezing and thawing dNTP;
5. Specific primers should be used and concentration of primers should be optimized and We suggest 0.4μM as a starting point for optimizing; Oligo-dT and Random primers are not suitable for this kit;

Required materials for RT

Instruments and consumables	Reagents
Centrifuge	DEPC
Pipettes	Tips
RNase free tubes	Electrophoresis Buffer
Water bath instruments	Loading Buffer
Gel electrophoresis instruments	DNA Marker

Reference


1. Houts, G.E., Miyagi, M., Ellis, C., Beard, D., and Beard, J.W. (1979) *J. Virol.* 29(2):517-522.
2. Guide to Molecular Cloning Techniques. Methods in Enzymology, Volume 152. pp 316-325. Edited by Shelby Berger and Alan R. Kimmel. Academic Press, Inc.

Description

BioRT cDNA First Strand Synthesis Kit provides the materials needed to rapidly and reliably synthesis first strand cDNA from RNA. It uses AMV reverse transcriptase to synthesize first strand cDNA from an RNA population.

The kit includes all the necessary reagents for first strand cDNA synthesis. The kit includes both random primers and oligo(dT18) primers. The user can choose either of these or alternatively use gene specific primers. The reverse transcriptase in the BioRT cDNA First Strand Synthesis Kit is AMV, which provides up to 60°C reverse transcriptase temperature and provides higher sensitivity and higher yield to cDNA's synthesis and PCR. The kit can synthesis up to 10kb cDNA. The reaction buffer is optimized for a kind of reaction system.

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