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MagaBio Blood Genomic DNA Purification Kit

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Appendix I

DNA purification evaluation

DNA yield is determined by measuring the concentration of DNA in the eluate by its absorbance at 260 nm. Absorbance readings at 260 nm should fall between 0.1 and 1.0 to be accurate.

Two formulations involved in DNA quantification is shown below:

Concentration of DNA sample = 50 µg/ml x A260 x dilution

Ratio=(A260-A320)/(A280-A320)

For accurate values, we recommend measuring absorbance in 10 mM Tris-HCl, pH 7.5.

Inhibitor and iron contaminant can be evaluated by the following PCR, real-time PCR, Southern, or other experiments.

Appendix II

The semi automation purification

With semi-automation machine, the Kit is deeply suitable for several samples, which supply a really platform of automation or streamline protocol and achieve high-throughput and high-speed but effective purification. An example for applying the kit on our product NPA-32:

▣ Sample processing

Add sample to 96 deep well plate. The processing method is the same to the above.

Note: When vortexing the plate please use special rotator, When Incubating at 56 °C please use special incubator.

▣ MagaBio adsorption

Add 250 µl of the Binding Buffer followed by 10µl of the **well-mixed** (particles are uniformly suspended) MagaBio Reagent. Vortex the plate for 30 seconds or so.

Note: when vortexing the plate please use special rotator, and the time need longer.

▣ Add 600µl washing buffer and 100µl elution buffer to the following well in turn, and add ddH₂O to the last well .then run the program.

Note: Please wash 2 or 3 times and elute 1 time in order to acquire the maximal result.

At one time the NPA-32 can process 32 samples at most, have many the default setting programs. If any interest please contact Bioflux Technical service immediately

Note: The OD value is probably higher than the anticipation because of the special method, which hardly impress the downstream application. You can estimate the result based on the electrophoresis. If any question, please contact Bioflux Technical service immediately.

Appendix III

FAQ

Q1: How does the MagaBio Nucleic Acid Isolation System differ from other existing magnetic products for the same applications?

A:

- I. A majority of magnetic-based nucleic acid isolation products use silica magnetic Particles / beads which they inherently have low capacity for nucleic acids.
- II. They appear to have low affinity for nucleic acids, since organic solvents such as ethanol are required during the washing steps.
- III. Some require two different wash buffers.
- IV. Other magnetic-based nucleic acid isolation products depend on antibody(ies) reaction or binding to other biological substance(s) that may require heat (up to 80°C) to inactivate the binders and liberate the nucleic acid. Also because of the requirement for an immunochemical reaction, they have limited applications.

Q2: What is the capacity of MagaBio magnetic particles?

A: 1 mg of MagaBio can bind up to 500 µg of calf thymus DNA with the Binding Buffer and without any chaotropic agents or detergents.

Q3: Does MagaBio Nucleic Acid purification System use an organic solvent?

A: No. The MagaBio Nucleic Acid purification kits reagents are environmental-friendly.

Q4: What is the expected typical DNA yield using MagaBio Blood Genomic DNA Purification kits?

A: Using MagaBio Blood Genomic DNA Purification Kit and 200 µl of fresh whole blood, the DNA yield of 5-12 µg (25-60 µg/ml) is expected.

Q5: What is the size of DNA isolated with the MagaBio DNA Isolation kit?

A: It depends on the application. With the MagaBio Blood Genomic DNA Purification Kit for genomic DNA isolation from i.e., blood, the purified DNA size is >20Kb.

Q6: Could DNA purified by MagaBio Blood Genomic DNA Purification Kit be save in frozen conditon?

A: Yes. After two years storage at -20°C, no significant loss of DNA has been observed.

Q7: Does MagaBio Nucleic Acid purification protocol require heating for elution of the purified nucleic acid?

A: No.

Q8: How many types of MagaBio DNA purification Kits are available? What applications are currently available?

A: There are only five types of MagaBio DNA purification kits that can be used to purification DNA from your entire sample source. The purification protocols for all three types of MagaBio kits are the same after the lysis step. This makes the MagaBio DNA purification Kit user-friendly and automation-friendly product.

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组成	数量
Protease K (PK)	1.0ml
Lysis Buffer	10.0ml
Binding Buffer	25.0ml
Wash Buffer	90.0ml*2
Elution Buffer	10.0ml
MagaBio Reagent	1.0ml
使用手册 V1.0	1 copy

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Important Notes

- Typically, >85% of the DNA is recovered in the first elution. If desired, more DNA can be recovered by applying a second elution.
- For RNA-free genomic DNA preparation, 2µl of an RNase A stock solution (10 mg/ml) should be added to the sample before addition of the Lysis Buffer.
- The DNA isolated from whole blood, using the conditions described in the Kit, typically shows an approximate size of 20-30 kb. DNA fragment of this length denatures completely during thermal cycling and can be amplified with the highest efficiency.
- Occasionally a tint of yellow color may be observed in the DNA isolated from whole blood. This will not affect the downstream processing such as PCR and sequencing.

Protocol

⊕ Sample processing

- Equilibrate all reagents and samples to room temperature.
- Pipet 10µl of PK Solution into the bottom of a 1.5 ml microcentrifuge tube.
- Add 100µl of sample to the microcentrifuge tube from the above.
- Add 100µl of the Lysis Buffer to the sample from the above and mix by pulse-vortexing for 15 seconds.
- Incubate at 56°C for 10 minutes.
- Remove the tube from 56°C.

⊕ MagaBio adsorption

- Add 250 µl of the Binding Buffer followed by 10µl of the **well-mixed** (particles are uniformly suspended) MagaBio Reagent.
- Mix the tube gently and incubate for 10 minutes at room temperature, while mixing.
Note: Using an end-over-end rotator or manual mixing every 2-3 minutes.
- Sediment the MagaBio DNA bound particles using a magnetic rack. Aspirate the supernate, remove the tube from the magnetic rack and wash the particles as described below.

⊕ Washing

- Add 500 µl of Wash Buffer to the tube from the above. Mix well by inverting the tube several times to ensure the particles are completely dispersed. Sediment the particles on the magnetic rack and aspirate the supernate.
- Remove the tube from the magnetic rack and repeat washing once more following the above step.

⊕ Elution

- Add 100~200µl of Elution Buffer and mix for 10 minutes.
Note: vortex gently every 2~3 minutes.
- Sediment the particles on the magnetic rack and carefully transfer the supernate containing the isolated DNA into a clean tube. The material is ready for further analysis. If the isolated DNA sample is not going to be tested on the same day, freeze at -20°C until the time of analysis.

Kit Components (100T)

Component	Amount
Protease K (PK)	1.0ml
Lysis Buffer	10.0 ml
Binding Buffer	25.0ml
Wash Buffer	90.0ml*2
Elution Buffer	10.0ml
MagaBio Reagent	1.0ml
Handbook V2.0	1 copy

Storage

- ▶ The Protease K is to be stored at 2~8°C, others at room temperature.
- ▶ All reagents, when stored properly, are stable for 12 months from the time of delivery.

Introduction

The Kit provides a very simple, fast and cost effective technique to isolate high quality DNA. Using one simple protocol, high yield of purified DNA can be isolated from various sources including whole blood, leukocytes cells. MagaBio sample processing is based on proprietary magnetizable particles--MagaBio. The pure DNA can be applied extensively in PCR, sequencing, Southern hybridization, mutant analysis, SNP and the others.

According to the special interaction, use MagaBio nucleic acid separation system with a general protocol---sample processing, MagaBio adsorption, washing and elution, and can go high-throughput.

Principle and Advantage

DNA in the sample is liberated using Protease K (PK) and Lysis Buffer. Released DNA is bound exclusively and specifically to the MagaBio Reagent in presence of Binding Buffer. The DNA bound to MagaBio particles is captured by the magnet and contaminants are removed by washing with Wash Buffer once or more. The DNA is then eluted from the particles with Elution Buffer or molecular grade water.

MagaBio Magnetic technical have great advantages:

- Mini sample, and high purification
- Simple and streamline separation procedure, used for auto-platform
- First elution can acquire 85% or more
- No organic solvent
- No high salt solution. no inhibitor
- No spin column, no centrifuge

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