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## MagaBio General Genomic DNA Purification Kit

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**Website: [www.bioer.com.cn](http://www.bioer.com.cn)**

sample is not going to be tested on the same day, freeze at -20°C until the time of analysis.

## Appendix I

### DNA purification evaluation

DNA yield is determined by measuring the concentration of DNA in the eluate by its absorbance at 260 nm. Absorbance readings at 260 nm should fall between 0.1 and 1.0 to be accurate.

Two formulations involved in DNA quantification is shown below:

Concentration of DNA sample = 50 µg/ml x A260 x dilution

Ratio=(A260-A320)/(A280-A320)

For accurate values, we recommend measuring absorbance in 10 mM Tris-HCl, pH 7.5.

Inhibitor and iron contaminant can be evaluated by the following PCR, real-time PCR, Southern or other experiments.

## Appendix II

### The semi automation purification

With semi-automation machine, the Kit is deeply suitable for several samples, which supply a really platform of automation or streamline protocol and achieve high-throughput and high-speed but effective purification. An example for applying the kit on our product NPA-32:

#### ▣ Sample processing

Add sample to 96 deep well plate. The processing method is the same to the above.

**Note: when vortexing the plate please use special rotator, When Incubating at 56°C please use special incubator.**

#### ▣ MagaBio adsorption

Add 250 µl of the Binding Buffer followed by 10µl of the **well-mixed** (particles are uniformly suspended) MagaBio Reagent. Vortex the plate for 30 seconds.

**Note: when vortexing the plate please use special rotator.**

#### ▣ Add 600µl washing buffer and 100µl elution buffer to the following well in turn, and add ddH2O to the last well .then run the program.

**Note: please wash 2 or 3 times and elute 1 time in order to acquire the maximal result .**

**At one time the NPA-32 can process 32 samples at most, have many the default setting programs. if any interest please contact Bioflux Technical service immediately**

Note: The OD value is probably higher than the anticipation because of the special method, which hardly impress the downstream application. You can estimate the result based on the electrophoresis. If any question, please contact Bioflux Technical service immediately.

## Appendix III

### FAQ

**Q1:** How does the **MagaBio** Nucleic Acid Isolation System differ from other existing magnetic products for the same applications?

**A: I** A majority of magnetic-based nucleic acid isolation products use silica magnetic Particles / beads which they inherently have low capacity for nucleic acids.

**II** They appear to have low affinity for nucleic acids, since organic solvents such as ethanol are required during the washing steps.

**III** Some require two different wash buffers.

**IV** Other magnetic-based nucleic acid isolation products depend on antibody(ies) reaction or binding to other biological substance(s) that may require heat (up to 80°C) to inactivate the binders and liberate the nucleic acid. Also because of the requirement for an immunochemical reaction, they have limited applications.

**Q2:** What is the capacity of **MagaBio** magnetic particles?

**A:** 1 mg of **MagaBio** can bind up to 500 µg of calf thymus **DNA** with the Binding Buffer and without any chaotropic agents or detergents.

**Q3:** Does **MagaBio** Nucleic Acid purification System use an organic solvent?

**A:** No. The **MagaBio** Nucleic Acid purification kits reagents are environmental-friendly.

**Q4:** What is the expected typical **DNA** yield using **MagaBio Genomic DNA Purification** kits?

**A:** i. Using **MagaBio General Genomic DNA Purification Kit** and 200 µl of fresh whole blood, the **DNA** yield of 5-12 µg (25-60 µg/ml) is expected.

ii. Using **MagaBio General Genomic DNA Purification Kit** purify various human tissues genomic **DNA**, **DNA** yields of 10-30 µg/5 mg lungs, spleen and cerebellum, and 25-65 µg/10 mg liver are expected.

**Q5:** What is the size of **DNA** isolated with the **MagaBio DNA** Isolation kit?

**A:** It depends on the application. With the **MagaBio General Genomic DNA Purification Kit** for genomic **DNA** isolation from i.e., blood, the purified **DNA** size is >20Kb. Plasmid **DNA** of smaller sizes (4~6 Kb) can be isolated with the **MagaBio General Genomic DNA Purification Kit**.

**Q6:** Could **DNA** purified by **MagaBio General Genomic DNA Purification Kit** be save in frozen condition?

**A:** Yes. After two years storage at -20°C, no significant loss of **DNA** has been observed.

**Q7:** Does **MagaBio** Nucleic Acid purification protocol require heating for elution of the purified nucleic acid?

**A:** No.

**Q8:** How many types of **MagaBio DNA** purification Kits are available? What applications are currently available?

**A:** There are only five types of **MagaBio DNA** purification kits that can be used to purification **DNA** from your entire sample source. The purification protocols for all three types of **MagaBio** kits are the same after the lysis step. This makes the **MagaBio DNA** purification Kit user-friendly and automation-friendly product.

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组成	数量
Protease K (PK)	0.5ml
RNase A stock solution	100µl
Lysis Buffer	5 ml
Binding Buffer	12.5ml
Wash Buffer	90ml
Elution Buffer	10ml
MagaBio Reagent	0.5ml
使用手册 V1.0	1 份

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## Important Notes

- Typically, >85% of the DNA is recovered in the first elution. If desired, more DNA can be recovered by applying a second elution.
- For RNA-free genomic DNA preparation, 2µl of an RNase A stock solution (10 mg/ml) should be added to the sample before addition of the Lysis Buffer.
- The DNA isolated using the conditions described in the Kit, typically shows an approximate size of 20-30 Kb. DNA fragment of this length denatures completely during thermal cycling and can be amplified with the highest efficiency.

## Protocol

### ☉ sample processing

- Equilibrate all reagents and samples to room temperature.
- Pipet 10µl of PK Solution into the bottom of a 1.5 ml microcentrifuge tube.
- Add 100µl of sample to the microcentrifuge tube from the above.  
**Note:**For cultured cells, use 200µl of appropriate number of cells (maximum 2.5 x 10<sup>7</sup>) resuspended in PBS.
- Add 100µl of the Lysis Buffer to the sample from the above and mix by pulse-vortexing for 15 seconds.
- Incubate at 56°C for 10 minutes.
- Remove the tube from 56°C.

### ☉ MagaBio adsorption

- Add 250 µl of the Binding Buffer followed by 10µl of the **well-mixed** (particles are uniformly suspended) MagaBio Reagent.
- Mix the tube gently and incubate for 10 minutes at room temperature, while mixing .  
**Note:** using an end-over-end rotator or manual mixing every 2-3 minutes.
- Sediment the MagaBio DNA bound particles using a magnetic rack. Aspirate the supernate, remove the tube from the magnetic rack and wash the particles as described below.

### ☉ washing

- Add 500 µl of Wash Buffer to the tube from the above. Mix well by inverting the tube several times to ensure the particles are completely dispersed. Sediment the particles on the magnetic rack and aspirate the supernate.
- Remove the tube from the magnetic rack and repeat washing once more following the above step.

### ☉ elution

- Add 100~200µl of Elution Buffer and mix for 10 minutes.  
**note:** vortex gently every 2-3 minutes.
- Sediment the particles on the magnetic rack and carefully transfer the supernate containing the isolated DNA into a clean tube. The material is ready for further analysis. If the isolated DNA

## Kit Components (50T)

Component	Amount
Protease K (PK)	0.5ml
RNase A stock solution	100µl
Lysis Buffer	5 ml
Binding Buffer	12.5ml
Wash Buffer	90ml
Elution Buffer	10ml
MagaBio Reagent	0.5ml
Handbook V1.0	1copy

## Storage

- ✎ The Protease K and RNase A is to be stored at 2~8°C, others at room temperature.
- ✎ All reagents, when stored properly, are stable for 12 months from the time of delivery.

## Introduction

The Kit provides a very simple, fast and cost effective technique to isolate high quality DNA. Using one simple protocol, high yield of purified DNA can be isolated from various sources including whole blood, buffy coat, leukocytes and cultured cells. MagaBio sample processing is based on proprietary magnetizable particles--MagaBio. The pure DNA can be applied extensively in PCR, sequencing, Southern hybridization, mutant analysis, SNP and the others.

According to the special interaction, use MagaBio nucleic acid separation system with a general protocol---sample processing, MagaBio adsorption, washing and elution, and can go high-throughput.

## Principle and advantage

DNA in the sample is liberated using Protease K (PK) and a Lysis Buffer. Released DNA is bound exclusively and specifically to the MagaBio Reagent in presence of a Binding Buffer. The DNA bound to MagaBio particles is captured by a magnet and contaminants are removed by washing with Wash Buffer once or more. The DNA is then eluted from the particles with an Elution Buffer or molecular grade water.

MagaBio Magnetic technical have great advantages:

- ☞ Mini sample, and high purification
- ☞ Simple and streamline separation procedure, used for auto-platform
- ☞ First elution can acquire 85% or more
- ☞ No organic solvent
- ☞ No high salt solution. no inhibitor
- ☞ No spin column, no centrifuge

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